



MARQMETRIX®

MEASUREMENT MADE SIMPLE

CONTACT US

BALLPROBE® - 0.25 IN.

0.25 in. (6.35mm) diameter, for small-diameter laboratory & light process applications with TouchRaman® immersion technology

DESIGNED FOR LAB & LIGHT-PROCESS SETTINGS

The MarqMetrix® 0.25 in. BallProbe is a smaller-diameter variant of our robust TouchRaman immersion contact probe.

Designed for moderate chemical environments in laboratory and light-process settings, this probe is constructed from stainless steel 316L, UV grade-sapphire, and an epoxy seal with wide chemical compatibility. The probe is designed to fit in cuvettes and other narrow access vessels, while providing the same optical performance as the Process BallProbe - 0.5 in. diameter.

SIMPLY TOUCH & MEASURE

MarqMetrix BallProbe technology utilizes an exclusively sourced, high-grade spherical sapphire lens. The short focal length of the spherical optic allows for TouchRaman—where users simply touch the probe to the sample—yielding highly reproducible sampling of liquids, solids, slurries, powders and heterogeneous mixtures.

The simplicity of design is especially important in process applications where measurement accuracy and reproducibility are mission critical.

**Reduced sampling variation,
extreme durability and
remarkable ease-of-use**



FEATURES & BENEFITS

- Immersion contact probe for touch measurements in moderate chemical environments
- Measures liquids, solids, slurries, powders, and heterogeneous mixtures
- Simple design for measurement accuracy and reproducibility
- Designed to fit cuvettes and other narrow access vessels

APPLICATIONS

- Biopharmaceutical manufacturing
- Pharmaceutical manufacturing
- Food and beverage processing
- Polymer and plastic manufacturing
- Hazardous chemical applications



For more information contact:

206-971-3625

info@marqmetrix.com

The power of Raman, the simplicity of MarqMetrix

Visit MarqMetrix.com



BALLPROBE - 0.25 IN.

0.25 in. (6.35mm) diameter, for small-diameter laboratory & light process applications with TouchRaman® immersion technology

Each probe accommodates a high-precision 6mm UV-grade sapphire ball lens and an internal focusing lens. The internal focusing lens improves performance and acts as an additional barrier against upstream chemical ingress.

The curvature of the sapphire ball facilitates material exchange near the surface of the lens, preventing the buildup of materials that interfere with spectral acquisition. The form factor and 'self-cleaning' properties make the BallProbe an ideal choice for process flow applications.

The utility of the BallProbe is optimized when paired with our filtered fiber-optic interface, creating the MarqMetrix Fiber BallProbe—a complete sampling solution for accurate and repeatable Raman measurements.



Specifications

Standard Probe Length	4 in. (102mm)
Probe OD (Outside Diameter)	0.25 in. (6.35mm)
Sample Working Distance	TouchRaman (Sample contacts BallProbe lens)
Continuous Operating Temperature Range	-20°C to 150°C
Pressured Design Condition	6,000psi (413 bar)
Compatible Wavelengths	500-1100nm

Operating Conditions

Suitable for intermittent exposure to most dilute and concentrated acids (hot & cold), bases and most organic solvents including ethanol, THF, ethyl acetate, acetone, DCM, toluene, pentane and acetonitrile

Avoid exposure to aqua regia, heated formic acid, mixed acid, TFAA

Wetted Materials

Probe Body	0.25 in. (6.35mm) OD Stainless steel 316L
Immersion Optics	6.00mm diameter UV-grade sapphire ball
Sealing Materials	High-temperature epoxy

Optical Properties

Made with high purity UV-grade sapphire ball lens aligned along the C-axis, eliminating response variability due to birefringence

Related Products

All-In-One Process Raman System - A single-unit spectrometer, laser, and acquisition computer

Process BallProbe® - 0.5 in. diameter, a larger diameter BallProbe for more demanding process or laboratory applications

Fiber BallProbe® - filtered fiber optic interface specifically designed for the MarqMetrix BallProbe and MarqMetrix Raman All-In-One

Proximal BallProbe® - Non-contact sampling probes for the analysis of in-motion samples through physical barriers

For more information call:

206-971-3625

info@marqmetrix.com

MarqMetrix.com

8016 v2.1

U.S. Patent No. 6,977,729 B2

